

ReSound Series: App Tips and Tricks

Enabling Bluetooth for your ReSound app

With the update to iOS 13 on your device, Apple® is now requiring users to consent to Bluetooth use in your apps outside of your phone's general Bluetooth.

If your ReSound app is asking for or requesting Bluetooth access, but your phone has Bluetooth turned on, tap **Settings** and scroll down and look for your ReSound app. Select the app and ensure that Bluetooth is enabled for your app.

	10:33	- I ? -	
	Settings Smart 3D		
	ALLOW SMART 3D TO ACCESS		
	Location While	Using >	
\langle	8 Bluetooth		>
	Microphone		
	Camera		
	💦 Siri & Search	>	
	Notifications Banners, Sounds, Badges	>	
	Background App Refresh		
	(1) Cellular Data		

